
Economical Magnetic Stirrer

Specification

Item	Parameter
Voltage [VAC]	100-240
Frequency [Hz]	50/60
Power Input [W]	1.2W
Adjustment	Scale
Speed Range [rpm]	300—2000
Max. Capacity[H ₂ O,L]	1.5
Max. Stir Bar	35mm
Work Platform[mm]	Ø120mm
Work Platform Material	PET
Casing Material	Anti-flaming-
Dimension D×W×H [mm]	145 ×150×50
Weight [kg]	0.3
Permissible Temperature [°C]	5-40
Permissible Humidity	80%
Protection Class	IP21

1. **POWER LED**
2. **SPEED ADJUSTMENT KNOB**
3. **WORK PLATFORM**
4. **SENSOR HOLDE**

Safety Precautions

The detailed precautions described in the User Manual are carefully developed in an attempt to cover all the possible risks. Please read carefully before operating this instrument.

- Ensure that only trained staff works with the instrument.
- Never use the instrument in any manner not specified in these instructions.
- If any abnormal condition occurs during operation, please stop it immediately and contact our service center.

Caution– Magnetism Effects of the magnetic field have to be taken into account (E.g. data storage media, cardiac pacemakers ...).

- Please pay attention to splash.

Operation

Placed the magnetic stirrer in the level of platform, confirm the magnetic stirrer is working under the rated voltage. After that, rotate the speed adjustment knob to control the magnetic stirrer stirring speed, clockwise is speed up, anti-clockwise is speed down. The magnetic stirrer can be stepless speed adjusting from 300rpm to 2000rpm.

Warranty


This product is guaranteed for two year from the date of delivery provided that they have been operated and maintained properly.

Troubleshooting and Error Code Information

Fault code	Problem description	Cause	Solution
E01	No operating response	No power supply	Check and connect the power supply before turning ON the power again.
		Power switch indicator lamp is OFF.	Switch power ON
		Connection failure	Open the instrument, check all connections and reconnect.
E02	Instrument does not stir	Target speed is not set	Set target speed, and ensure that the speed indicator is ON
		Driver board failure	Replace the driver board
		Motor failure	Replace the motor component

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