
Mini Dry Bath

Heating Mini Dry Bath

Heating & Cooling Mini Dry Bath

Heating & Cooling Mini Dry Bath with Hot Lid

Please carefully read the Instructions, and use the product safely under the direction of the Instructions.

Disclaimer

- The intellectual property right and relevant copyright belongs to the company.
- Products of the company are under the patent protection of China and other countries and regions (including obtained patents and patents being applied for).
- We elaborately prepared the manual in the attitude of being responsible for the users. But we can't guarantee that the contents of the manual are fully correct. If any occasional or subsequent loss is caused by the use of the Manual, the company will not be liable for this at all.
- The manual is a pure technical document, free of any implication or oblique hint of any third party. Moreover, we won't be liable for any user's misunderstanding upon the printer's error(s).
- The company and any of the company's employees will not bear any liability for direct/indirect loss of information or stop of business (if any) caused by the information of the manual or the product it mentions.
- The company reserves the right to change specification and price of the product.
- If any information of the manual changes, we will not notify this otherwise.
- The manual and any content of it must not be copied, abstracted or modified in any form without the company's written approval in advance.

Table of Contents

1. SAFETY NOTES	1
2. Product overview.....	2
2.1 SCOPE OF APPLICATION	2
2.2 TECHNICAL PARAMETERS	2
2.3 Components	3
3. OPERATING INSTRUCTIONS	4
3.1 OPEN-PACKAGE INSPECTION.....	4
3.2 OPERATION.....	4
4. CLEANING AND MAINTENANCE.....	8
5. FAULT DIAGNOSIS	9
PRODUCT CERTIFICATION.....	9
WARRANTY.....	10

1. Safety notes

 **Caution!**

- Please carefully read the Manual before operating the product, and observe the specifications on safe operation.
- The product can be operated only by one who has been specially trained.

 **Caution Hot!**

- Please touch the base and the heating module with care. The maximum temperature of the heating module is as high as 100°C (specific temperature shall be subject to technical parameter).
- After the shutdown of the product, care should be taken not to be scalded by the afterheat on the heating module.

 **Ground protection for safety!**


- For the purpose of safety, please confirm the power socket has been properly grounded before putting the product into service.

- Please make sure that the power voltage meets the voltage requirement specified on the nameplate.
- Please use the product indoors. Don't use it in dangerous environment or in water.
- The instrument shall be put at a place with low humidity and less dust, far away from water source and without the perpendicular incidence of sunshine and intense light source. The room shall be well-ventilated, free of corrosive gases or strong magnetic interference, and far away from heating gas, fire and all other heat sources.
- There are air vent holes at the bottom and both sides of the device. During the operation of device, these air vents must be not blocked or shielded by any object but well-ventilated.
- Before each startup of the product, please confirm that the product and its accessories are undamaged. For the purpose of safety, please use the standard accessories listed in the chapter "Accessories" in accordance with the Manual. Accessories must be firmly connected to the product in a way that avoids separations.
- During the work of the product, please wear proper protective equipment. You need to make clear the possibility of and be held responsible for the dangerous circumstances above by yourself.
- The power cable must be kept away from the surface of heating module rather than shield the device.
- Please prevent water from splashing on the electrical elements of the product.
- The power supply can be fully disconnected only if the power plug is pulled out. Before assembly, disassembly, cleaning or maintenance, the main power plug must be pulled out.

2. Product overview

2.1 Scope of application

- The product is intended for heating liquids in the environments like school, laboratory and factory. An environment where the product works should meet the following requirements.
- Voltage fluctuation: within the range of -10%~+10% of normal value (the product is designed for indoor socket);
- Min. distance between instruments, instrument and wall:

The product is unusable in residential area or under the restrictions specified in Chapter 1.

Fig. 2.1

2.2 Technical parameters

Model	Heating & Cooling Mini Dry Bath with Hot Lid	Heating & Cooling Mini Dry Bath	Heating Mini Dry Bath
Display	LCD	LCD	LCD
Range of temperature setting [°C]	0~100	0~100	25~100
Range of temperature control [°C]	Room temperature -23~100	Room temperature -23~100	Room temperature +5~100
Accuracy of temperature control [°C]	±0.5	±0.5	±0.5
Accuracy of temperature adjustment [°C]	0.1	0.1	0.1
Heating time (25°C~100°C)	≤20min	≤20min	≤20min
Cooling time (20°C~0°C)	≤25min	≤25min	/
Maximum heating rate	8°C /min	8°C /min	6.5°C /min
Maximum refrigeration rate	3°C /min	3°C /min	/

Range of timing	0~999min/ 0~999sec	0~999min/ 0~999sec	0~999min/ 0~999sec
Program	9 pcs (2 stages for each)	9 pcs (2 stages for each)	9 pcs (2 stages for each)
Fast calibration	Supported	Supported	Supported
USB interface protocol	Supported	Supported	Supported
Error code hint	Supported	Supported	Supported
External dimension [L x W x H mm]	110x162x140	110x162x14 0	110x162x140
Weight of whole machine [kg]	≤1	≤1	≤1
Power supply [v]	100~240V, 50/60Hz	100~240V, 50/60Hz	100~240V, 50/60Hz
Maximum Power [w]	60	60	60
Operating temperature [°C]	+10~40	+10~40	+10~40
Operating humidity[% rH]	≤80	≤80	≤80

2.3 Components



Hot cover (for Heating & Cooling Mini

Dry Bath with Hot Lid type only)

Heating module

Display screen

Operation
panel

Power switch Power line

USB interface

Fig. 2.3.1

Fig. 2.3.2

Name	Description
ON/OFF	Start / stop program
Prog.	Select a program

▲ ▼	Set time / temperature parameter
➡	Enter parameter setting / move the cursor

Table 2

3. Operating instructions

3.1 Open-package inspection

Open-package inspection aims to confirm the completeness of associated parts. The packing list is given below:

Name	Qty.
Main unit	1
Power line	1
Power adapter	1
Handle	1
Operating manual	1

If you find any damage on the package, please specify the damage on the receipt. If you find any internal damage after opening the package, please contact local supplier or manufacturer.

Caution:

If you find any visible damage on the product, please don't connect the product to power supply.

3.2 Operation

3.2.1 Preparation

- The device shall be steadily placed on a work platform, especially a ventilated place without flammables or explosives.
- Inspect whether the grid voltage falls within the working voltage range specified on the nameplate;
- Make sure the power socket is well-grounded;
- Connect power line well, and open power switch.

3.2.2 Create a program

- To press the “Prog.” key, and select a program;
- For example, select P1, and then press “➡” key to enter the temperature/time parameter setting interface, later press “➡” to move the cursor and “▲” or “▼” key to modify temperature/time value.
- Every program may set two stages, “■” means temperature and time of the first stage, and “■” means temperature and time of the second stage.

P1	37	5min	■
P1	60	3min	■

Temperature and time of the first stage

Temperature and time of the second stage

The second stage may be closed and not used. For setting the time of the second stage, the “▼” key shall be pressed for a while till “OFF” appears, by then the second stage will be closed.

P1	60	OFF	■
----	----	-----	---


3.2.3 Modify the time unit

Please keep pressing the “Prog.” key and turn on the power switch, and please not release the “Prog.”Key until the setting interface of time unit pops up, and then press the “➡” key to select time unit between minute (min) and second (sec), finally press the “ON/OFF” key for confirmation.


3.2.4 Run a program

- To press the “Prog.” Key, and select one of the 9 Preset Programs P1,P2...P9.
- For example, after selecting Preset Program P1, the “ON/OFF” key may be pressed to run it, and the device will start to run the first stage of program. When the temperature rises, “NOTOK” will be displayed on the right side of display screen; when the temperature reaches the set value, a long beep will sound to prompt this, and the countdown will begin, and the set time, e.g. “3 min”, will be displayed on the right side of display screen, and the icon “■” of the first stage of program will flash.


If the time of the first stage of program is set as “0min”, when the temperature of the first stage reaches the set temperature, a long beep will sound to prompt this, the

constant temperature state will begin, and “Hold” will be displayed on the right side of display screen.

- After the timekeeping of the first stage ends, the device will automatically run the second stage of program, similarly, when the temperature rises, “NOTOK” will be displayed on the right side of display screen; when the temperature reaches the set value, a long beep will sound to prompt this, and the countdown will begin, and the set time, e.g. “1 min”, will be displayed on the right side of display screen, and the icon “■” of the second stage of program will flash.

After the timekeeping of the second stage of program ends, a long beep will sound to prompt this, and the display interface will turn back to the first stage.

3.2.5 Stop a program

During program running, if the “ON/OFF” key is pressed, the program will stop and the timekeeping will end.

Normal work of the product during the operations above means the product can be put into normal service. If abnormal work occurs, it means the product may be damaged during transport, please contact the after-sales service center of local supplier or manufacturer.

3.2.6 Restore the factory settings

Please keep pressing the “➡” key and turn on the power switch, and please not release the “➡” key until the “System Reset?” interface pops up in the display screen, and press the “ON/OFF” key to confirm the restoration of factory settings.

3.2.7 Temperature calibration

Before delivery, the temperature of device has been calibrated, but the actual temperature may deviate from the real temperature for some reasons, and the temperature deviation may be corrected by temperature calibration.

It is recommended that the temperature must be calibrated every time based on temperature calibration function after change of module, in order to ensure temperature control accuracy.

The device adopts the three-point temperature calibration method, i.e. synchronous linear calibration at the three temperature points 37°C, 60°C, 90°C. The linear calibration of the three temperature points guarantees that the accuracy of the three temperature points of the system within $\pm 0.5^{\circ}\text{C}$.

At the time of temperature calibration, the ambient temperature must be 15 to 25 °C.

Caution:

Before delivery, the temperature of device has been calibrated, so please don't use this function at will!

Calibration procedure:

- To fill silicone oil in the center hole of module, and put the probe of thermometer in the hole (The thermometer shall be accurate to 0.1°C, and its probe must be fully immersed in the hole.);
- To keep pressing the “ON/OFF” key and turn on the power switch, please not release the “ON/OFF” key until the “System Cali?” interface pops up in the display screen, finally press the “ON/OFF” key to confirm the entry to the calibration interface;
- The device will first automatically control temperature at 37°C, and remain at the constant temperature for 20min, then the reading of thermometer shall be read, if the reading is 36.5°C, the “▲” or “▼” key shall be pressed to modify the displayed value in the screen to be 36.5, finally the “ON/OFF” key shall be pressed for confirmation;

P0	37	36.2
----	----	------

- After the confirmation by the “ON/OFF” key, the device will automatically rise to 60°C, and remain at the constant temperature for 20min, then the reading of thermometer shall be read, if the reading is 60.5°C, the “▲” or “▼” key shall be pressed to modify the displayed value in the screen to be 60.5, finally the “ON/OFF” key shall be pressed for confirmation;

P0	60	60.3
----	----	------

- After the confirmation by the “ON/OFF” key, the device will automatically rise to 90°C, and remain at the constant temperature for 20min, then the reading of thermometer shall be read, if the reading is 90.5°C, the “▲” or “▼” key shall be

pressed to modify the displayed value in the screen to be 90.5, finally the “ON/OFF” key shall be pressed for confirmation;

P0	90	90.2
----	----	------

- After the calibration of all the three points, if the display screen alternatively shows the “Cali Successful” interface and the K, B value interface, it can be considered that the calibration ends; if the display screen only shows the “Cali Failed” interface, the input value of calibration may be wrong, and the temperature calibration shall be repeated once in this case.

4. Cleaning and Maintenance

Operate and maintain the product properly, so that it is in a good working state, which can extend the service life of the product. In routine service, keep the product dry and clean, remove the spilled liquid quickly, clean the outer surface with a non-grinding cleaner, and do not connect the power supply until all surfaces are dry. If liquid or moist solid enters the product, please disconnect the power supply quickly and leave off, and contact the manufacturer / supplier for more advice.

- The deep hole of module shall be regularly cleaned by clean soft rag with a little of anhydrous alcohol so that the test tube fully contacts the hole wall to realize good heat conduction and avoid pollution.
- Stains on the surface of device, if any, shall be cleaned by clean soft rag with detergent.
- Keep the product clean, and the cleaning solution is not allowed to flow into the machine.
- Power must be disconnected before maintenance and cleaning, and please use our recommended methods to clean the product.
- It’s forbidden to clean the module or the device by corrosive cleaning solution.
- When the product is not used for a long time, please store the product with power off and place it in a dry, clean and smooth place in normal temperature.

Warning!

Before any maintenance or inspection, the power cable must be pulled out of the socket.

Please refer to *Service Manual* for details on machine maintenance.

5. Fault diagnosis

The advanced production technology and testing method are used for the dry bath and a rigorous testing is conducted for each instrument before delivery, so that this product has good reliability. In service, the common failures are generally caused by improper operation or setting. If there are faults that cannot be handled, please record the fault phenomenon and notify the local dealer, or contact us directly.

The common phenomena of the dry bath in the application are shown as follows:

Fault phenomena	Causes	Solutions
Error 1	Fault of circuit board or sensor	Return to depot for repair
Error 2	Fault of hot cover circuit board or sensor	Return to depot for repair
Error 3	The displayed temperature exceeds the temperature control range.	1. Restore factory settings; 2. Return to depot for repair
Error 4	Hot cover temperature exceeds the temperature control range.	1. Restore factory settings; 2. Return to depot for repair
The display screen does not turn bright after power-on	Loose connection of power cable	To inspect whether the power cable is firmly connected
	Poor contact of switch	Change switch
	Controller damage	Return to depot for repair
The module does not get hot or cold	Damage of internal circuit	Return to depot for repair

Product certification

Our company has certified that the product meets applicable national product standards and industrial product standards of China as well as ISO9001 standards, and further certified that the product meets applicable standards published by other members of the ISO.

Warranty

According to the manufacturer's warranty clause, the warranty period of the product is 24 months provided that it is used by the normal operation method under the normal service conditions specified in the Manual. If a problem covered by the manufacturer's warranty clause occurs, please contact local supplier. Or, you can directly mail the product to the manufacturer with the packing list and a problem description attached. You need to bear the transportation expense incurred by such mailing.

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